



Raymond Price^{1*}, Pauravi J. Gandhi¹, Christian Petroski¹, Dong Liu¹, Chen Tian¹,
Vincent Truong², Patrick Walsh², Robert Petroski¹

¹Neuroservices-Alliance, San Diego CA 92121, ²Anatomic Inc., Minneapolis, MN 55413




Abstract-




- Human iPSC-derived nociceptors are an increasingly important model system for drug discovery assays in pain. However, it is important to characterize the cell properties to understand and qualify iPSCs as a 'fit-for-purpose' model.
- We recorded from over 100 Anatomic RealDRG™ hiPSC-derived nociceptors from 1-8 weeks in vitro. Human nociceptors plated at low density on rat astrocytes appeared healthier and were easier to patch compared to neurons plated without astrocytes. Nociceptors grew larger with time in culture as evidenced by an increase in Cm, with decreased Rm. The RMP became more hyperpolarized after three weeks and remained stable thereafter.
- We observed little spontaneous activity in cells at every developmental timepoint, as would be expected from sensory neurons. The nociceptors were excitable, however, as action potentials were elicited with depolarizing current injections. Further, the amount of current injection required to elicit the first action potential (rheobase) increased over developmental time, as expected from neurons with lower Rm. Additionally, nociceptors at all developmental timepoints exhibited sodium currents.
- The voltage-dependent sodium current amplitudes increased with time in culture, while the peak sodium current density (pA/pF) remained stable. We also evaluated the presence of TTX-R sodium currents at 5-8 weeks in vitro and detected their presence at all evaluated time points. The proportion of TTX-R among nociceptors varied. We also observed an increased spontaneous neuronal activity in sensory neurons with longer time in vitro, using HD-MEA.
- Overall, we determined that Anatomic RealDRGs® develop a mature neuronal phenotype by three weeks in culture, with hyperpolarized resting membrane potential, the ability to evoke spike trains, and exhibit TTX-sensitive and resistant Na currents, suggesting they may be suitable for further assay development.

Methods-

- Anatomic RealDRGs® were thawed and plated at a low density, with a subset plated on a monolayer of rat astrocytes.
- Whole cell patch clamp recordings were performed on human iPSC-derived nociceptors at different timepoints (1-8 weeks) to determine the time course for the development of a mature neuronal phenotype.



- Upon establishing the whole cell configuration, we measured the passive membrane properties for every cell. We also recorded spontaneous action potential firing, intrinsic excitability as well as evaluated the expression of voltage-gated sodium currents and tetrodotoxin (TTX)-resistant sodium channels (in presence of 0.5 μM TTX).
- We measured spontaneous action potentials using the Maxwell MaxTwo high-density multielectrode array (HD-MEA) instrument.



Conclusions-

- Human iPSC sensory neurons exhibited a mature developmental phenotype by 3 weeks *in vitro*, hyperpolarized resting membrane potential and ability to evoke action potentials.
- Human iPSC sensory neurons elicited an increase in spontaneous neuronal firing as they mature in culture.
- These sensory neurons expressed key proteins and transcripts.
- These sensory neurons showed both TTX-S and TTX-R Na currents. Na current amplitude after TTX was large enough (~1000 pA) to test Nav 1.8 blocking compounds.

